

Infectin Sample Protocols

Overview – Infectin™ is a viral infection enhancer designed to facilitate viral penetration of the cortical actin barrier, thereby greatly enhancing productive viral infection. Infectin™ can be used to facilitate the infection of a variety of host cells by different viruses and viral vectors. Below is an example procedure for enhancing lentiviral infection of suspension CEM-SS T Cells with Infectin™. We recommend using this procedure as a starting point and testing multiple doses to determine the optimal concentration for your cells.

Suspension Cell Protocol – Lentiviral infection enhancement of suspension CEM-SS T-cells

- 1) Count cells and determine viability.
Note: Cell viability should be $\geq 80\%$.
- 2) Pellet cells by centrifugation at 300 x g for 5 minutes.
- 3) Resuspend cells in complete media at a concentration of $\sim 2 \times 10^6$ cells/mL.
- 4) Add 100 μ L of cell suspension ($\sim 2 \times 10^5$ cells) per infection.
- 5) Pre-treat cells by adding 10 μ L of Infectin™ (10X) so that Infectin™ concentration is approximately 1X.
- 6) Mix and incubate for 10-30 minutes at 37°C in a CO² incubator.
- 7) Thaw viral particles on ice.
- 8) Add the desired volume of virus for infection to a separate microcentrifuge tube.
- 9) Add the volume of Infectin™ (10X) equal to 10% of the volume of the virus to the microcentrifuge tube containing the virus.
- 10) Mix gently.
- 11) Add the virus/Infectin™ mixture to the cells and incubate the infection for 2 hours at 37°C in a CO² incubator.
- 12) Add 1 mL of fresh complete media.
- 13) Centrifugation at 300 x g for 5 minutes.
- 14) Remove supernatant
- 15) Add 1 mL of fresh complete media.
- 16) Culture the infected cells for 2-3 days.
- 17) Quantify viral infection.

Table 1: Scaleup recommendations for viral infection using Infectin™

Infection Cell Number	Cell Volume	Infectin™ Buffer (10X)	Final Volume
2×10^5	100 μ L	10 μ L	1 mL
5×10^5	250 μ L	25 μ L	2.5 mL
1×10^6	500 μ L	50 μ L	5 mL
2×10^6	1 mL	100 μ L	10 mL
5×10^6	2.5 mL	250 μ L	25 mL
1×10^7	5 mL	500 μ L	50 mL
5×10^7	25 mL	2.5 mL	250 mL
1×10^8	50 mL	5 mL	500 mL

Infectin Sample Protocols

Adherent Cell Protocol – Lentiviral transduction enhancement of adherent HDFn cells

- 1) Count cells and determine viability.
Note: Cell viability should be $\geq 80\%$.
- 2) Pellet cells by centrifugation at 300 x g for 5 minutes.
- 3) Resuspend cells in complete medium at a concentration of $\sim 2 \times 10^6$ cells/mL.
- 4) Add 250 μ L of cell suspension ($\sim 5 \times 10^5$ cells) to a 6-well plate.
- 5) Add 1750 μ L of complete medium to each well.
- 6) Gently mix and incubate for 4-12 hours at 37°C in a CO² incubator.
- 7) Remove medium.
- 8) Wash cells with 2mL of fresh medium.
- 9) Remove 1.5mL of wash medium (0.5mL medium remaining).
- 10) Pre-treat cells by adding 50 μ L of Infectin™ (10X) so that Infectin™ concentration is approximately 1X.
- 11) Mix and incubate for 10-30 minutes at 37°C in a CO² incubator.
- 12) Thaw viral particles on ice.
- 13) Add the desired volume of virus for infection to a separate microcentrifuge tube.
- 14) Add the volume of Infectin™ (10X) equal to 10% of the volume of the virus to the microcentrifuge tube containing the virus.
- 15) Mix gently.
- 16) Add the virus/Infectin™ mixture to the cells.
- 17) Mix and incubate the infection for 2 hours at 37°C in a CO² incubator.
- 18) Remove medium
- 19) Add 2 mL of fresh complete media to wash.
- 20) Remove wash medium
- 21) Add 2 mL of fresh complete media.
- 22) Culture the infected cells for 2-3 days.
- 23) Quantify viral infection.