

Infectin™ Sample Protocol

Overview – Infectin™ is a viral infection enhancer designed to facilitate viral penetration of the cortical actin barrier, thereby greatly enhancing productive viral infection. Infectin™ can be used to facilitate the infection of a variety of host cells by different viruses and viral vectors. Below is an example procedure for enhancing lentiviral infection of suspension CEM-SS T Cells with Infectin™. We recommend using this procedure as a starting point and testing multiple doses to determine the optimal concentration for your cells.

(If Infectin™ is received in powder form, please add 700µl of the provided Infectin™ Dissolving Buffer. Mix gently to avoid generating bubbles. Once resuspended, Infectin™ can be stored at 4°C for 3 months)

Suspension Cell Protocol – Lentiviral infection enhancement of suspension CEM-SS T-cells

- 1) Count cells and determine viability.
Note: Cell viability should be ≥ 80%.
- 2) Pellet cells by centrifugation at 300 x g for 5 minutes.
- 3) Resuspend cells in complete media at a concentration of ~2 x 10⁶ cells/mL.
- 4) Add 100µL of cell suspension (~2 x 10⁵ cells) per infection.
- 5) Pre-treat cells by adding 10 µL of Infectin™ (10X) so that Infectin™ concentration is approximately 1X.
- 6) Mix and incubate for 10-30 minutes at 37°C in a CO² incubator.
- 7) Thaw viral particles on ice.
- 8) Add the desired volume of virus to cells.
- 9) Add additional Infectin™ (10X) equal to 10% of the volume of virus added to the cells to maintain an Infectin™ concentration of 1X.
- 10) Mix gently.
- 11) Infect for 2 hours at 37°C in a CO² incubator.
- 12) Add 1 mL of fresh complete media.
- 13) Centrifugation at 300 x g for 5 minutes.
- 14) Remove supernatant.
- 15) Add 1 mL of fresh complete media.
- 16) Culture the infected cells for 2-3 days.
- 17) Quantify viral infection.

Table 1: Scaleup recommendations for viral infection using Infectin™

| Infection Cell Number | Cell Volume | Infectin™ Buffer (10X) | Final Volume |
|-----------------------|-------------|------------------------|--------------|
| 2 x 10 ⁵ | 100 µL | 10 µL | 1 mL |
| 5 x 10 ⁵ | 250 µL | 25 µL | 2.5 mL |
| 1 x 10 ⁶ | 500 µL | 50 µL | 5 mL |
| 2 x 10 ⁶ | 1 mL | 100 µL | 10 mL |
| 5 x 10 ⁶ | 2.5 mL | 250 µL | 25 mL |
| 1 x 10 ⁷ | 5 mL | 500 µL | 50 mL |
| 5 x 10 ⁷ | 25 mL | 2.5 mL | 250 mL |
| 1 x 10 ⁸ | 50 mL | 5 mL | 500 mL |

Adherent Cell Protocol – Lentiviral transduction enhancement of adherent HDFn cells

- 1) Count cells and determine viability.
Note: Cell viability should be $\geq 80\%$.
- 2) Pellet cells by centrifugation at 300 x g for 5 minutes.
- 3) Resuspend cells in a complete medium at a concentration of $\sim 2 \times 10^6$ cells/mL.
- 4) Add 250 μ L of cell suspension ($\sim 5 \times 10^5$ cells) to a 6-well plate.
- 5) Add 1750 μ L of complete medium to each well.
- 6) Gently mix and incubate for 4-12 hours at 37°C in a CO² incubator.
- 7) Remove medium.
- 8) Wash cells with 2mL of fresh medium.
- 9) Remove 1.5mL of wash medium (0.5mL medium remaining).
- 10) Pre-treat cells by adding 50 μ L of Infectin™ (10X) so that Infectin™ concentration is approximately 1X.
- 11) Mix and incubate for 10-30 minutes at 37°C in a CO² incubator.
- 12) Thaw viral particles on ice.
- 13) Add the desired volume of virus to cells.
- 14) Add additional Infectin™ (10X) equal to 10% of the volume of virus added to the cells to maintain an Infectin™ concentration of 1X.
- 15) Mix gently.
- 16) Mix and incubate the infection for 2 hours at 37°C in a CO² incubator.
- 17) Remove medium
- 18) Add 2 mL of fresh complete media to wash.
- 19) Remove wash medium.
- 20) Add 2 mL of fresh complete media.
- 21) Culture the infected cells for 2-3 days.
- 22) Quantify viral infection.