

# Rapid Cell-Based SARS-CoV-2 Neutralization Assay with HA-CoV-2 Pseudovirus

## Firefly Luciferase Reporter and READY TO USE CELLS

**Purpose:** Neutralization Assay utilizing the HA-CoV-2 pseudovirus for titering antibody containing serum using multiple dilutions and 18hr incubation on a 96 well plate. **Equipment. Supplies and Reagents** 

Item	Source, Catalog #, Lot #
Samples	•
Human antibody (AB) serum	Patients or purchased AB solutions
Equipment Not Included	
Sterile Reagent Reservoir	Fisher #03-449-255
Olympus 96-Well PCR Plate, Non-Skirted Ultra Thin Wall, Natural, 25 Plates/Unit	Genesee Scientific # 24-300
96-Well Tissue Culture Plate, Greiner Bio-One (With Lid, μClear® White Flat Round, Chimney)	VWR #82050-758
Luciferase Reader: GloMax® Discover Microplate Reader	Promega #GM3000
Reagents Included	•
HEK293T(ACE2+TMPRSS2)#21 cells READY TO USE FROM THE CRYOTUBE TUBE	Virongy Biosciences Inc. # RCSNAK-01
VBI Infection Medium	Virongy Biosciences Inc. # RCSNAK-01/ RCSNAK-02
HA-CoV-2(Luc) pseudovirus	Virongy Biosciences Inc. # RCSNAK-01/RCSNAK-02
Standard SARS-CoV-2 Neutralizing Antibody	Virongy Biosciences Inc. # RCSNAK-01/RCSNAK-02
10X Cell Lysis Buffer	Virongy Biosciences Inc. # RCSNAK-01/RCSNAK-02
D-luciferin Solution	Virongy Biosciences Inc. # RCSNAK-01/RCSNAK-02
Firefly Luciferase Buffer Solution	Virongy Biosciences Inc. # RCSNAK-01/RCSNAK-02



VIRONGY BIOSCIENCES INC.

### Firefly Luciferase Reporter and READY TO USE CELLS

### **Neutralization Assay Procedure:**

\*All cell culture work should be performed in a laminar airflow Class II biosafety cabinet and strict sterile procedures must be used.

#### 1. Prepare the Plate Maps: Mixing Plate (96-well PCR Plate):

Recommended Mixing Plate Map						
	Dilution 1	Dilution 2	Dilution 3	Dilution 4	Dilution 5	Dilution 6
	1:5	1:25	1:125	1:625	1:3125	1:15625
Sample 1	10µL AB	10µL Dilution 1	10µL Dilution 2	10µL Dilution 3	10µL Dilution 4	10µL Dilution 5
_	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium
Sample 2	10μL AB+40μL	10µL Dilution 1	10µL Dilution 2	10µL Dilution 3	10µL Dilution 4	10µL Dilution 5
	Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium
Sample 3	10µL AB+40µL	10µL Dilution 1	10µL Dilution 2	10µL Dilution 3	10µL Dilution 4	10µL Dilution 5
_	Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium
Sample 4	10μL AB+40μL	10µL Dilution 1	10µL Dilution 2	10µL Dilution 3	10µL Dilution 4	10µL Dilution 5
_	Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium
Sample 5	10µL AB+40µL	10µL Dilution 1	10µL Dilution 2	10µL Dilution 3	10µL Dilution 4	10µL Dilution 5
	Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium
Sample 6	10µL AB+40µL	10µL Dilution 1	10µL Dilution 2	10µL Dilution 3	10µL Dilution 4	10µL Dilution 5
	Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium
Sample 7	10µL AB+40µL	10µL Dilution 1	10µL Dilution 2	10µL Dilution 3	10µL Dilution 4	10µL Dilution 5
	Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium	+40µL Medium
Standard	10  uI A $B + 40  uI$	10uI Dilution 1	10uI Dilution 2	10uI Dilution 3	10ul Dilution 4	10uI Dilution 5
Neutralizing	Medium	+40ul Medium	+40ul Medium	+40ul Medium	+40ul Medium	+40ul Medium
Antibody	wiedłum			+opt Medium	- τομε Medium	

## Infection Plate (96-well White Chimney TC Treated Plate):

Recommended Infection Plate Map													
	1	2	3	4	5	6		7	8	9	10	11	12
Sample 1	1:20	1:100	1:500	1:2500	1:12500	1:62500	Sample 1	1:20	1:100	1:500	1:2500	1:12500	1:62500
Sample 2	1:20	1:100	1:500	1:2500	1:12500	1:62500	Sample 2	1:20	1:100	1:500	1:2500	1:12500	1:62500
Sample 3	1:20	1:100	1:500	1:2500	1:12500	1:62500	Sample 3	1:20	1:100	1:500	1:2500	1:12500	1:62500
Sample 4	1:20	1:100	1:500	1:2500	1:12500	1:62500	Sample 4	1:20	1:100	1:500	1:2500	1:12500	1:62500
Sample 5	1:20	1:100	1:500	1:2500	1:12500	1:62500	Sample 5	1:20	1:100	1:500	1:2500	1:12500	1:62500
Sample 6	1:20	1:100	1:500	1:2500	1:12500	1:62500	Sample 6	1:20	1:100	1:500	1:2500	1:12500	1:62500
Sample 7	1:20	1:100	1:500	1:2500	1:12500	1:62500	Sample 7	1:20	1:100	1:500	1:2500	1:12500	1:62500
Standard							Controls	IC	IC	IC	PBC	PBC	PBC
Neutralizing	1:20	1:100	1:500	1:2500	1:12500	1:62500							
Antibody													

Infections: 15µL AB + 45µL HaCoV2(Luc) Particles + 15µL cells

IC= Infection Control (15µL Medium + 45µL HaCoV2(Luc) Particles + 15µL cells)

PBC= Particle Background Control (30µL Medium + 45µL HaCoV2(Luc) Particles)

#### 2. Antibody Treatment Preparation in a 96-Well PCR Plate:

\*\*The volumes are doubled in the mixing plate for each sample to distribute on the **Infection Plate** in duplicate. In a 96-well PCR plate (**Mixing Plate**), prepare the antibody(AB) serum dilutions according to the Mixing Plate map.

- 1. Place 40µL of VBI Infection Medium in each row of the Mixing Plate.
- Perform a serial dilution with each of the serum sample and control standard antibody according to the table below: Add 10 μL of each AB/serum sample to the first column and mix well. Transfer 10 μL from column 1 to column 2 and mix well using a higher volume setting on the pipette, pipette-mix the AB/serum and medium thoroughly before transferring to the next



## Rapid Cell-Based SARS-CoV-2 Neutralization Assay with HA-CoV-2 Pseudovirus

### Firefly Luciferase Reporter and READY TO USE CELLS

column. Continue to dilute until the 6<sup>th</sup> dilution is complete. **Note:** Be sure to exchange pipette tips between well transfers and ensure that antibody/serum and medium are mixed thoroughly to produce accurate results.

Well #	Antibody Serial Dilutions
1A	10 μL of AB/Serum +40μL VBI
(1:5) AB	Infection Medium
2A	$10 \ \mu L \text{ of well } #1 + 40 \ \mu L \text{ VBI Infection}$
(1:25) AB	Medium
<b>3A</b>	$10 \ \mu L \text{ of well } #2 + 40 \ \mu L \text{ VBI Infection}$
(1:125) AB	Medium
4A	$10 \ \mu L \text{ of well } #3 + 40 \ \mu L \text{ VBI Infection}$
(1:625) AB	Medium
5A	10 $\mu$ L of well #4 + 40 $\mu$ L VBI Infection
(1:3125) AB	Medium
6A	$10 \ \mu L \text{ of well } \#5 + 40 \ \mu L \text{ VBI Infection}$
(1:15625) AB	Medium

- 3. In the **Infection Plate add 45 μL** of HA-CoV-2(Luc) particles to each infection well following the Infection Plate map.
- 4. Transfer 15  $\mu$ L of the AB/serum dilutions from the **Mixing Plate** to the **Infection Plate** with the HACoV2 particles and mix well. Repeat this step to prepare duplicate wells according to the Infection Plate map.
- 5. Set 3 wells as the **Infection Control** wells by adding **15µL** of VBI Infection Medium to the particles and mix well
- 6. Set 3 wells as the **Particle Background Control** wells by adding **15μL** of the VBI Infection Medium to the particles and mix well. **Note: These wells will not receive cells in the next step.**
- 7. Incubate the HA-CoV-2(Luc) particles and AB mixture for 30 mins at 37°C. Prepare the cells while the particle and antibody incubate.

<b>Final Dilutions</b> Antibody/Serum + Virus Complex					
Α	1:20				
В	1 : 100				
С	1 : 500				
D	1 : 2500				
E	1 : 12500				



Rapid Cell-Based SARS-CoV-2 Neutralization Assay with HA-CoV-2 Pseudovirus

VIRONGY BIOSCIENCES INC.

Firefly Luciferase Reporter and READY TO USE CELLS



# 1. Prepare the HEK293T(ACE2+TMPRSS2) READY TO USE CELLS

- 1. Ouickly thaw the HEK293T(ACE2+TMPRSS2)#21 READY TO USE cells in a 37°C water bath. Note: water bath can quickly become contaminated, its recommended to wrap the tube with parafilm.
- 2. Once thawed immediately transfer the tubes contents to 5 mL of prewarmed VBI Infection Medium.
- 3. Centrifuge the cell for 5 minutes at 1200rpm and discard the supernatant. Resuspend the cells in 2mL of prewarmed VBI Infection Medium.
- 4. Transfer 15µL of the resuspend cells into each infection well of the Infection Plate. Each well should receive around 2.5x10<sup>4</sup> cells in 15 µL of VBI Infection Medium. Note: For the Particle Background Control wells add 15µL of VBI Infection Medium.
- 5. Incubate the Infection Plate overnight at 37°C for 16-18 hours. Note: This incubation can be done in the GloMax plate reader and can be set to automatically add the lysis buffer when completed with the injector add on.

# 2. Luciferase Procedure

- 1. Following the incubation of the Infection Plate add 7.5 µL of the 10X Cell Lysis Buffer directly to each well and mix by orbital shaking for 2 mins. Allow cells to lyse for at least 5 minutes.
- 2. In a dark location, away from direct light, prepare the Firefly Luciferase Assay Solution by mixing the D-Luciferin Solution with the Firefly Luciferase Buffer Solution in a 1:50 ratio. For a full 96 well plate prepare 5.5 mL of the luciferase substrate solution by mixing 110 µL of the D-Luciferin Solution with 5.390 mL of the Substrate Buffer. Note: Use the Firefly Luciferase Assay Solution within 30 minutes of preparation. Do not reuse the Firefly Luciferase Assay Solution.
- 3. Add 50 µL of the Firefly Luciferase Assay Solution to the cell lysates. Note: For more accurate reading the Firefly Luciferase Assay Solution should be injected into each well by a programable injector so the timing of the assay can be precisely controlled.
- 4. The plate should be run within 10-15 minutes of adding the samples to the plate. (Tip: Run the plate reader in a dark room to reduce any background signal.) Note: More lysate can be added to generate a stronger luciferase signal. Calculate the results by subtracting the background signal determined from the average of the "Control No Infection" wells.